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Radiation Sterilization Verification Guide for Disposable Bioprocessing Systems

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Disposable biological process manufacturing systems are increasingly being implemented in the biopharmaceutical industry based on the benefits of safety, time saving and cost reduction. These disposable systems are used to process or contain liquids from culture media, additives and buffers, to bulk semi-finished products and final formulations. In many cases, control of microorganisms and sterility is required to ensure the purity and safety of the product. Radiation sterilization is a common microbiological control and sterilization method suitable for use in disposable systems. In the pharmaceutical industry, the standard method of verifying radiation sterilization is not widely known, and the industry has historically relied on temperature and humidity (steam) to sterilize stainless steel bioprocess systems.

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Under the auspices of the Plastics Industry Association, the BPSA Technical Committee established a subcommittee to develop guidelines for irradiation and sterilization of disposable biological process components and systems.

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Scope of application: The purpose of this guide is to provide basic information and recommendations for the use of gamma radiation for microbial control and sterilization in the biopharmaceutical industry.

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Although electron beam radiation is similar to gamma ray radiation and follows the same criteria, it is generally not used for the microbiological control or sterilization of single-use systems due to its limited penetrating power. Because of its limited penetrating ability, single-use systems, therefore It is not considered in this document.

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Uses: The BPSA Radiation Sterilization Subcommittee proposes the following in this guide:

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biopharmaceutical industry, make relevant experiments that are sufficient or required for applications in different biopharmaceutical production processes Education decision.

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• Summarize the standard method of sterilization of medical products with gamma radiation, which is recognized by the biopharmaceutical industry for verification of single-use systems.

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• Provides recommendations for bioburden and sterility testing methods used in validation of radiation sterilization of biological process components and systems.

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Product radiation has established organizations such as the American National Standards Institute (ANSI), which promotes medical devices (AAMI), the International Organization for Standardization

(ISO) and ASTM International. Advances in the biopharmaceutical and biotech drug industries have led to the need for pre-sterilized or microbially controlled products that can be directly introduced in key manufacturing processes. Initially, these disposable products were targeted for small-scale applications such as laboratory-scale drug development and preclinical studies. In the last few years, single-use disposable product lines at the clinical and process scales, including high zone filter capsules, membrane chromatography units, tubing, connectors, and flexible biological vessels (eg, polymer membrane bags) With an exponential increase. For radiation sterilization, verifying that established industry-standard medical products have been successfully applied to such products, sterility verification needs and equipment and drug regulatory agencies recognize it. The sterilization validation standard for medical product biotechnology applications, however, can be quite expensive and complicated. Because of the further discussion in this guide, the use of the master mode of the system, through the standard support, can reduce these burdens, or no radiation Sterility testing is used where possible, where a high but undefined microbial control is sufficient. Key definitions are presented in the "Glossary".

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Basic knowledge of gamma radiation

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Gamma radiation is the application of electromagnetic radiation (rays) emitted from radioisotopes such as cobalt 60 (^{60}Co) and cesium 137 (^{137}Cs) isotopes. Gamma rays are not obstructed by most materials and can penetrate most disposable bioprocess system components. Ionizing radiation causes the inactivation of microbial nucleic acids by destroying nucleic acids. Gamma rays are also not retained by matter, leaving no residual radioactivity. The unit of gamma radiation dose is kGy, which is the energy absorbed by its quantified radiation. A gray scale is the amount of 1 Joules of radiation energy absorbed by one thousand gallons of material (one thousand Gore = 1 Joule/Kelvin). A dose of ≥ 8 kGy is usually sufficient to eliminate low bioburden levels (Reference 2 in Table 9). In an elevated bioburden level (>1000 colony forming units or cfu per unit) higher doses may be required to achieve sterility. Typically, a sterile 10^{-6} Aseptic Assurance Level

(SAL) is achievable at 25 kGy. Even with higher bioburden levels, bioburden can be reduced with a lower probability of achieving sterility (for example, SAL 10^{-5} or 10^{-4}). Irradiate this SAL

The level of product is still sterile, but the probability of non-sterility is still high, may not meet the requirements of the industry standard to confirm the requirements of the sterilized standard sterilized medical products gamma irradiation process using well-defined operating parameters to ensure accurate dose. In a well-designed irradiation facility, the only variable for any given density of material determines the material cost in the radiation field when the irradiated product and the amount accepted by the microorganism are determined. The product is not exposed to heat, humidity, pressure or vacuum. The minimal byproducts of waste generated by gamma irradiation and the need for isolated degassing (such as sterilization with ethylene oxide gas) or routine biological reaction tests. As a constant and predictable sterilization method, gamma irradiation is safe, time-saving and cost-effective. In addition to microbial inactivation, gamma radiation also causes ionizing polymer molecules and excitation (7). Over time and depending on the dose, the absorbed dose can result in changing the physical or chemical properties of the polymer. Some polymers exhibit higher resistance than other radiation-induced changes, but all polymers are affected. While the minimum dose of 25 kGy required for microbial control and sterilization, the actual dose applied is usually in the 25-50 kG range. It is to be noted that the adsorption of radiation is cumulative and that some polymers may be suitable for use in gamma irradiation. The use of one-off biological process systems may not be able to withstand more than 50 kGy and still maintain their integrity over time. Therefore, we should avoid repeated irradiation (8).

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Microbial control and sterilization

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The general concept of medical product sterilization protocols described by current gamma irradiation in current industry standards can be directly applied to disposable biological process systems. These industry standards require the validation of the effects and reproducibility of the sterilization process, based on which the average bioburden and the system's minimum radiation dose are irradiated after the subsequent sterility test. The validated sterile system is also subject to routine audits involving bioburden and sterility tests. Standard procedures were performed for gamma sterilization validation and applied separately to the 6th and 7th disposable biological process systems. These standards for pre-clinical or clinical drug development applications in small-scale custom

batch manufacturing, particularly for use in a sample system biological treatment system, it can be expensive and cumbersome, and final product design sterilization re-verification may be required. As an alternative to aseptic labeling requirements, many biological process systems can be controlled by appropriate microorganisms, typically at 25 KGy doses under irradiation, and sterilization verification by standard methods may not be necessary. The minimum 25-kilo-Geie dose sterilization claim originated from Charles Artandi and Wall. The authors of a study conducted by Dun Van Winger in 1959 identified "lowest dose killing" of over 150 different microorganisms. As a result of their research, they chose 25 kGy for the sterilization dose, stating that "[25 kGy] is the least resistant to killing more than 40% of the most intractable microorganisms" (9). As a result, 25 kGy was established as a suitable minimum exposure for sterilization. The historical sterilant dose of 25 k Gei can also be sufficient to eliminate surviving bioburden and provide high levels of microbial control without the need for validated sterility requirements. Components that require zero or low bioburden when applied in a non-sterile process do not require verification of sterility requirements and may be a qualified system as a microbial control. The following sections are intended to help readers determine whether to verify sterilization, or more simply, microbial control, is appropriate

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Bio-load control in the use of single-use system bio-load control stand-alone

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In biopharmaceutical production, where the drug product is produced by cell culture, it can be divided into four stages of treatment (Figure 1). In the upstream phase, nutrients and other fluids, with the production of cells, are introduced into a fermenter or bioreactor and the cells produce the target molecule. During the harvest phase, the cells are separated by centrifugation of target molecules, deep filtration or membrane filtration. In the downstream stage, the target molecule undergoes a series of separations, purifications, and chromatography and membrane filtration to produce the final, purified raw material that reaches the drug product concentration level. In the final formulation and filling phase, the final liquid formulation in the purified whole (active pharmaceutical ingredient (API) or biological) is prepared, filter sterilized, or filled into a sterile final container under aseptic conditions. The latter stage is similar to how synthetic medicines are manufactured using a sterile process. In the course of these phases, there is a continuous need for microbial control to prevent undesirable pollution. International regulations state that aseptic processing must be validated to support sterilization requirements. This does not apply. To be treated as microbiological control, it is also possible to have a zero or low bioburden, but it has not yet been

demonstrated that there is a defined SAL. Although sterility is guaranteed to pass the verification to provide the required raw materials as final sterilized preparations and to inject injected medicines and biological products, the preparation and disinfection of many manufacturing stages in biological processes (Fig. 1) is difficult to sterilize and may not require actual verification. Sterile. Instead, they are provided according to microbiologically controlled operations, such as verification by gamma irradiation without sterilization. Consider whether a process step must be validated as sterile, with a key contrast in process development between the controlled microbiological (zero or low bioburden) states and the time and cost requirements for optimization. Similarly, when verification of sterility is required, consideration of whether a complete disposable system is required or only the fluid path is confirmed as sterile can also affect the time and cost requirements. As mentioned earlier, establishing and maintaining sterility as a single use of the validation requirements for components or components of a biological process system can require that their users ultimately bear significant costs. Especially in the subsequent process of process development and microbiological control, the irradiation requirements can be sufficient to prevent significant microbial contamination, while saving the cost of maintaining an effective sterile requirement. The balance between increased costs to achieve this assurance provides a balance between supplemental protection is a must for biopharmaceutical manufacturers. When the choice must be made, it is most likely to be driven by the risk assessment of each unit operation and the manufacturing process as a whole. In the upstream processing, bacterial (e.g., *E. coli*) cell cultures tend to operate in a short period of time, and are low levels of bacteria that over-grow well tolerant contaminants. In this case, a user may think that a one-time use of a single use system to prepare disposable bioreactors and culture media or other nutritive additives can be safely handled according to the irradiation microbe

control provided without sterilization verification. Mammalian cell cultures, by contrast, are generally more susceptible to bacterial contamination than extended periods of time (eg, two run weeks). In such cases, verification is provided by sterilization to provide additional assurance that effective testing, although aseptic control is enhanced by γ -irradiation, the uncertain probability can be developed and even enhanced during the production process. Especially in the course of development cell harvesting and downstream processing steps (centrifugation, depth or membrane filtration, and chromatography, with intermediate holdings) are rare. Instead, the process equipment is chemically disinfected or sterilized and maintained as a microorganism controls zero or low bioburden. The middle shelf can be kept at a low temperature to prevent microbial growth and a non-sterile treatment step. In the event that these unprotected process steps are sterile, it is not necessary to verify the preparation of the sterile disposable system buffered feed solution or intermediate containment. Buffered solutions or intermediates are filtered through an irradiating bio-load reduction filter system. Microbial growth control is generally applied to process steps that are not validated as sterile. For the formulation and filling also offers the option to choose between radiation-sterile validation and control of microorganisms. Subsequent processing, including formulation, sterilization filtration, and ultimately aseptic filling of the aseptically finished active pharmaceutical ingredient in bulk or at a reduced temperature of biological storage, may be sufficient in a microbially controlled container prior to pre-filter storage. The only key step that clearly stipulates a single-use system is identified as sterility in preparation of sterile API and aseptically filled sterile containers; that is, if the fluid is required to be sterile, a disposable system It is populated to have verification requirements.

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Current standards for sterilization verification

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Some industry standards can be used for sterilization verification of gamma radiation health products. The ANSI / AAMI / ISO 11137:2006 standard (sterilized health supplements - radiation) was originally published in 1994, followed by a variety of supporting technical reports developed by ISO and AAMI. The key in these are the AAMI Technical Information Report (TIR) 33: Substantialization of selected sterilization doses in 2005 - Method VD Max (10) In ISO System Review 11137:1994 (US/ANSI/AAMI/ISO passed in the United States 11137:1994), the file was modified and the overall heading is divided into three sections, Sterilization of Health Products - Radiation's Three Parts of ANSI/AAMI/ISO 11137: Descriptive Than 2006. Part 1: Requirements for routine control of medical devices during development, validation, and sterilization processes dictate the need for radiation-sterilized healthcare products for development, verification, process

control, and routine monitoring. Part 1 applies to continuous and batch gamma irradiators using radioisotope ^{60}Co or ^{137}Cs , and uses a beam to irradiate electrons or X-ray generators. Part 2: Determination of Sterilization Dose The described method can be used to determine the minimum dose required to achieve the specified requirements for sterility, including methods to confirm 15 or 25 kGy as the sterilization dose. Part 3: Guided dosimetry aspects provide guidance doses for health care products, and establish maximum dose (product certification) of dosimetric aspects of radiation sterilization; establishment of sterilization doses; installation qualifications; business qualifications; The ANSI/AAMI/ISO standard describes two methods for sterilization validation:

- 剂量设定的方法 (方法 1 或方法 2, ANSI / AAMI / ISO 11137), 它考虑到分布和产品生物负载的抗辐射性, 并且在有限的程度, 最终使用的产品, 以及
 - The dose setting method (Method 1 or Method 2, ANSI/AAMI/ISO 11137), which takes into account the radiation resistance of the distribution and product bioburden, and to a limited extent, the final use of the product, and
 - 剂量属实方法, (VD 最大值的方法), 这意味着实验设计来限定预定的伽马剂量作为灭菌剂量 (1-3)。
 - Dosage is a faithful method, (method of maximum VD), which means that the experiment is designed to define a predetermined gamma dose as the sterilization dose (1-3).
- 其中一些方法'相似之处包括:
- Some of these methods' similarities include:
- 只有结合适当的质量管理体系的应用
 - Only apply the application of an appropriate quality management system
 - 承认的剂量控制和生物负荷控制是两个同等重要的方面保证无菌
 - Admitted dose control and bioburden control are two equally important aspects of aseptic assurance
 - 要求通过定期测试论证持续有效性以验证建立辐照灭菌剂量保持适宜产品灭菌。
 - It is required to demonstrate continuous

effectiveness through regular tests to verify the establishment of irradiation doses to maintain proper product sterilization.

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Table 1 summarizes the additional similarities and differences in the various methods. Although this illustrates different main attribute methods, please refer to the actual ANSI/AAMI/ISO standards or AAMI TIR for comprehensive information. In general, these methods provide the user with a degree of application to define the flexibility of the irradiation process standard to indicate whether the product is manufactured in all cases and commercialized under uniform conditions. Users are given scope to make decisions Be realistic, as long as they adhere to the correctness of the technology in several key areas:

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• Frequency of audit sampling (periodically to support the continued use of established sterilization doses). Although the initial development of the AAMI TIR33:2005 and ANSI/AAMI/ISO 11137 standards focused on the medical device industry, they also apply to systems and components of biological processes. The rest of this guide is therefore based on these industry standards.

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The dose setting method requires experimentation that results in the production of a sterile product with a qualified minimum radiation dose. 2006 and include: The following steps are described in ANSI/AAMI/ISO11137 when using the general method of this type of method:

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Get the product specified number of tests from the three production batches (Table 1).

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Experiment to determine the average bioburden under ANSI/AAMI/ISO11737-1 (11-12).

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Perform dose verification studies (irradiation and sterility tests).

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septic processing: Technically, disease-free microorganisms. In bioprocessing applications, it refers to a controlled environment in which an operation is performed or a connector designed with an AND gate to prevent contamination by introduction of microorganisms (which can be sterilized or microbiologically controlled).

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Bioburden: Active microbial population in a product. In the case of sterilization, it should be immediately determined before disinfection. Bioburden can describe its entire content in a project, including the outer surface, or only the fluid contact surface.

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Bioprocessing components: single-use components, modules, or a system that can include, but are not limited to, filters, membrane chromatography units, tubing, connections, connectors, flexible bags, or rigid containers and probes.

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Bioprocessing system: A combinatorial component of a combinatorial or biological process designed to facilitate the handling of drugs or single-use disposable formats in biological solutions. A bioburden reduction/microbiology system can be pre-assembled with control steps or assembled by the user by connecting pre-sterilized by using a sterile connection assembly.

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Dose Release: This product is a sterile assay based on physical irradiation processing data, not sterility

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False Positive Aseptic Test Results: Test results are incorrectly interpreted as indicating that the test article has a viable organism and the article is not sterile. Typically, such results result from the failure of the treated sample or from the interaction of the medium with the test or from the interaction between a sample and the growth of microorganisms.

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Irradiation: The application of ionizing radiation can destroy the bacteria content.

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Microbial control: The product is exposed to gamma radiation in a controlled clean assembly environment. This process reduces the bioburden load but does not support "sterile" labeling requirements.

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Sample Part (SIP): A defined part or component of a product assembly.

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Sterile: Does not contain live microorganisms.

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Sterility Assurance Level (SAL): Aseptic Assurance Level (SAL): In practice, sterility cannot be demonstrated. In contrast, sterility is expressed as a probability that one single item of live microorganisms will sterilize (1). Usually this term is expressed as 10^{-n} , not as an absolute value.

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Bacteria test: Test to determine if viable microorganisms are present. Typically, the test involves impregnating the assembly or system or rinsing the fluid path with a non-sterile microbial growth medium, the lower incubation conditions of said medium facilitating the growth of microorganisms, and observing a suitable other indication latency of turbidity or microbial growth.

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Sterilization: Used to render products from live microorganisms using proven methods. It is generally believed that the terminal disinfection unit appears to be aseptically achieving a level of aseptic assurance $\leq 10^{-6}$; that is, less than or equal to one chance at a probability of one million, and a viable microorganism is present in the sterilization unit (1-6). Lower SALs can verify aseptic (2) in some cases.

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Verification: Establishing documented evidence that provides a high level of assurance that a particular process will continue to produce products that meet their intended specifications and quality attributes.

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Figure 1: Schematic representation of a typical biological process using virtually all disposable equipment (exceptions include centrifuges, columns,

and flow filters or depth filters for tangential stainless steel holders); upstream and formulation, and fill in the equipment used in the step may require sterility verification (indicated by dashed lines), but the equipment has a more general harvest and the downstream stage may need to be irradiated only at the ≥ 25 kGy dose required for microbial control, and sterilization verification (Shown in solid or dashed lines).

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Interpret results - Accept or reject the study.

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Establish minimum dose disinfection based on the bioburden and SAL requirements.

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Regularly study the effectiveness of the audit dose test. The dose setting method involves the dose separating the product based on the actual radiation resistance of the qualifying microorganism. As a result, these methods provide the lowest possible dose, but require more product samples than the dose-based method.

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The ANSI/AAMI/ISO 11137-2:2006 VD

Maximum Methodology allows for the sterilization of two pre-established irradiation sterilization doses: 15 kGy and 25 kGy. AAMI TIR33: 2005 provided the flexibility to use seven additional sterilization doses: 17.5, 20, 22.5, 27.5, 30, 32.5, 35 kGy. Sterilization doses range from 15 kGy bioburden up to 1.5 colonies, 25 kGy bioburden up to 1,000 CFU, and 35 kGy bioburden up to 440,000 cfu. Manufacturers should maintain strict bioburden control at each stage of their production process. In addition to the supplemental dose setpoints, ANSI/AAMI/ISO 11137-2:2006 and AAMITIR33:2005 are based on the same principles and provide basically the same guidance. The general method when using the VD max method includes the following steps:

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法 1 和 2 的 ANSI / AAMI / ISO 11137 以下。Regularly study the effectiveness of audit dose tests. In the VD method, the largest demonstration of irradiation sterilization dose is 10^{-6} SAL; this method is not supported by the current higher SAL (lower guaranteed sterility). Two inputs for fixed-dose keying off: bioburden (experimentally determined) and worst case bioload resistance data driving assumptions (2). Although the choice of AAMI TIR33 Sterilization Dose: 2005 is fixed, sampling requirements are less than methods 1 and 2 of ANSI/AAMI/ISO 11137.

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Implementation method verification

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Given their size and complexity, verification tests in the process of one-off biological process system and component construction are difficult. A large-scale system can be, for example, a pipe connected to several meters, represented by a high zone filter capsule (s), which can be a 10,000-liter bag assembly of 254 mm (10 inches) or more, not a complicated connector. Design - or the complex manifold of the connector. Therefore, you should carefully plan a simple and scientifically effective method for bioburden and sterility testing. Sample Item Division (SIP) product groupings, and general methods in the ANSI/AAMI/ISO11137 test description: The 2006 guidelines provide strategies to address the design of medical devices and the complexity of medical products, and also apply to large-scale biological process components and such

as The system described in the sections described below.

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The simplest strategy is to sample three standard production lots. Typically, one or more product families developed by multiple product manufacturers not a single product may represent the status of the rest of their system (10). Used to indicate a product system. This item may be designated as the "main product," an "equivalent product," or "simulation product."

Main product: The main product is not necessarily the largest one. Where all materials and components are the same, and the product changes only in the size or quantity of the material, the main product is usually the largest one. Depending on the materials contained and the manufacturing process, the main product can be a combination of manufacturing and assembly, or a combination of these factors, the most typical operation of the assembly, the largest combination of materials, and the highest filtration area. In the end, the main product is a sterilization process that is considered the most challenging.

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The equivalent product: The equivalent products are not necessarily those who have the same bioburden level because the bioburden level

varies between the natural sample and the splice tube time sample. Whether the product can be considered equivalent;

- 生物负载估计是差不多的;
- The bioburden estimate is similar;
- 生物负载特性是相似的;
- Bioburden characteristics are similar;
- 生物负载估计结果在灭菌剂量是大约相同的使用方法 1 和方法 2，或验证剂量是大约相同的使用 VDmax 时。
- The bioburden estimate results at approximately the same sterilization dose as Method 1 and Method 2, or verify that the dose is approximately the same when using VDmax.

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Analog Products: A simulated product is a component that contains the same material and actual product but has been manufactured for testing purposes only and is therefore not an actual product. This is usually the worst-case property that is made to include or combine product lines. The composition or structure with the simulated product does not affect the correctness of the dose determination. In designing a simulated product, care should be taken to create a laboratory worker who does not produce undue challenges and who performs bioburden and sterility tests.

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Large-scale systems for implementing pathway testing

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Bioburden and bioprocess systems or there are two specific challenge components for sterility testing. First, large products are particularly difficult to handle sterile. Specific issues include the introduction of such large items into biosafety cabinets, the addition and removal of fluids to test items under aseptic conditions, the large volume of test fluid required to operate under sterile conditions, and the provision of incubation space for large systems, The additional problems that are used to test the entire product (rather than a separate fluid path) and treat its exterior surfaces without contaminating them must be overcome. The second problem is that biological process systems or components are often quite expensive. Taking such a product is sold at a very low volume at the beginning of its life cycle, and testing a relatively large number of actual items can be prohibitively commercialized.

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The specific method that manufacturers take to solve these problems must be developed specifically for each case to maintain a balance between them

- 一个愿望， 并以确保技术的正确性通过测试实际的文章，
- A desire, and to ensure the correctness of the technology by testing the actual article,
- 为了避免错误结果的愿望（虚增生物负载估计或假阳性无菌试验结果）所造成的不实用操作制品。
- Improper handling of products caused by the desire to avoid erroneous results (inflated bioburden estimation or false positive sterility test results).

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Several common methods can alleviate the difficulties caused by these two problems. The following examples are intended to illustrate what problems can be overcome; they are not intended to exclude other methods of scientific defense.

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Sample Project Section (SIP): This SIP method allows the manufacture and testing of a scaled product. This reduces the cost of test articles and the difficulty of manipulating products in the lab. This method adds experiments to determine bioburden estimates by multiplying the size between SIP-based articles and the actual product represented. The SIP multiplier can determine the largest system manufactured by experiment or simply by taking the ratio of the model.

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Fluid Path: In some cases, it is appropriate to verify the sterility of the product's fluid path. Testing the fluid path can be better than testing the entire product because one product is easier

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20 International Biopharmaceutical Engineering Process May 2008 Supplement

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A fluid path generally requires partial filling of the product with a sterile buffered bioburden test to

ensure that all surfaces are wetted, and then the product is agitated to manually push the suspended organism into the buffer zone. The buffer was then removed and the load of the bioburden by standard microbiological methods tested. To better quantify bioburdens, bioburden recovery studies were conducted before actual bioburden studies were initiated. The recycling study has two purposes: to quantify the efficiency of the bio-load recovery from the test strip, and to help determine the specific method that will be used in the bio-load study.

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For sterility testing, the article is either partially filled with a microbial growth medium or with an extract. In the following stirring process, the growth medium may be left inside or removed, and the microbial load may be measured.

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Segmented large-scale products: In order to facilitate the operation, a large number of products can be divided into components that are easier to

manipulate in laboratory tests or irradiation. This method is particularly effective if the entire product, rather than simply the fluid path is tested. An example of such a method is to cut a large biological container into 10 parts, cut into multiple small pieces each and add a piece of the container, a microbiology laboratory and a biological safety cabinet that can easily be handled under aseptic conditions. In this kind of In case, each container will be tested and the sum of the bioburden will represent the original bag's bioburden. Similarly, all 10 parts need to be sterile as a raw biological container counted as sterile.

The product packaging is used as a containment device: it is generally packaged in two bioprocess systems or components to ensure its cleanliness and sterility. Some products, especially those that have to be tested in their entire product, are very suitable for testing while they remain in the original packaging. In this case, the inner package serves as a storage container. A small opening is included in the package entering the test article and adding or removing a buffer or growth medium. The openings are more sterile than sealed to maintain a closed system and prevent microbial contamination. This method is less effective. The test article is a closed loop system that must be opened to access the interior surface.

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In the microbiological control or sterilization cycle of the gamma radiation, there must be means of quantifying the dose of gamma radiation. This is done by dosimetry, which quantifies the field's absorption of radiation. The measuring device is called a dosimeter. There are two aspects of dosimetry that are particularly important:

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sterilization dose must be accurately measured ($\pm 10\%$)

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- The load should map the dose to ensure that the uniformity of the gamma dose is sufficient to configure the radiation.

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Although dose behind science is beyond the scope of this guide, various dosage techniques exist. The user should ensure that the dosimeter and the operator are qualified to provide sufficient accuracy and accuracy, especially within the range of radiation dose considered. This method validates the concept of application, as it will be the data on which products are certified and released and can be yielded according to any measurement method. Typically, gamma-irradiated specialized service providers perform qualification dosimeters.

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After the irradiation cycle is verified, the microbes in production can control the use of sterile products. The product was released based on dosimetry after gamma irradiation: a dosimetric verification of the cycle produced the required minimum dose of product load. No routine sterility tests require significant release. Due to the product release dose, its manufacturer must ensure that the doses administered by the gamma radiography service providers have good conditions and that the gamma irradiation operation controls consistently

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Disposable biological systems represent a significant benefit to the biopharmaceutical process development and manufacturing. The benefits are accompanied by the adoption of these systems and the integration of new decisions by all stakeholders. With respect to quality attributes such as microbiological control and sterilization, some problems may be solved in general, while others should be handled in a specific biopharmaceutical production process. Figure 1 provides decision tree sterilization and microbial control.

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cThis guideline discusses the basic issues related to the level of manufacture of systems and components for the implementation of one-time biological processes for gamma irradiation and microbial control. It is designed to frame the problem of microbial control and make sensible, rational decisions to provide the structure at the level of microbiological control required for the manufacture and use of disposable biological process systems.

REFERENCES [参考文献]

- 1 ANSI/AAMI/ISO 11137-1:2006. Sterilization of Health Care Products — Radiation, Part 1: Requirements for Development, Validation and Routine Control of a Sterilization Process for Medical Devices; www.iso.org.
- 2 ANSI/AAMI/ISO 11137-2:2006. Sterilization of

Health Care Products — Radiation,
Part 2: Establishing the Sterilization Dose;
www.iso.org.
3 ANSI/AAMI/ISO 11137-3:2006. Sterilization of
Health Care Products — Radiation,
Part 3: Guidance on Dosimetric Aspects;
www.iso.org.
4 DIN EN 556. Sterilization of Medical Devices —
Requirements for Medical Devices to Be
Designated —Sterile; <http://webstore.ansi.org>.
5 US Pharmacopeia (USP). <1211>, Sterilization
and Sterility Assurance;
www.usp.org; www.uspnf.com/uspnf/login.
6 ANSI/AAMI ST67:2003. Sterilization of Health
Care Products — Requirements for Products
Labeled —Sterile; www.ansi.org.
7 Hemmerich, K.J., Polymer Materials Selection for
Radiation-Sterilized Products.MDDI, February
2000;
www.devicelink.com/mddi/archive/00/02/006.html.
8 AAMI TIR17. Radiation Sterilization — Material
Qualification;
www.aami.org; <http://webstore.ansi.org/ansidocstore/product.asp?sku=AAMI+TIR17%3A1997>.
9 Tallentire, A. The Spectrum of Microbial
Radiation Sensitivity. Radiat. Phys. Chem.15, 1980:
83–89.
10 AAMI TIR33:2005. Sterilization of Health Care
Products — Radiation,
Substantiation of a Selected Sterilization Dose –
Method VDmax;
www.aami.org; <http://webstore.ansi.org/ansidocstore/product.asp?sku=AAMI+TIR33%3A2005>.
11 ANSI/AAMI/ISO 11737-1:2006. Sterilization of
Medical Devices — Microbiological Methods, Part
1: Determination of a Population of
Microorganisms on Products; www.iso.org.
12 ANSI/AAMI/ISO 11737-2:1998. Sterilization of
Medical Devices — Microbiological Methods, Part
2: Tests of Sterility Performed in the Validation of a
Sterilization Process (first edition);
(下面是作者和对本文写作提供帮助和支持的人)
The members of the Irradiation and Sterilization
Subcommittee, formed under the umbrella of BPSA,

are corresponding author **Suraj Baloda**
(subcommittee chair), corporate quality manager at
Millipore Corporation, 290 Concord Road, Billerica,
MA 01821; 1-781-533-2641, fax 1-781-533-3395,
suraj_baloda@millipore.com; www.millipore.com;
Jerold Martin (chair, technology committee),
senior vice president of Pall Life Sciences; **Jeff
Carter**, director, filtration R&D at GE Healthcare;
Ernie Jenness, senior engineer at Millipore
Corporation; **Brett Judd**, validation engineer at
Thermo Fisher Scientific; **Kent Smeltz**, president
of Advanced Scientifics Inc.; **Isabelle Uettwiller**,
regulatory affairs director for Sartorius Stedim
Biotech SA; and **Melissa Hockstad**, senior
technical director for the Society of the Plastics
Industry, Inc.